




**Product Specification**  
**Aspiration pipette, unpadded, 2 ml, sterile,  
non-pyrogenic/endotoxin-free, non-cytotoxic, 1  
piece(s)/blister**

**Page 1**



**Product description**

Order number	86.1252.011
Product description	Aspiration pipette, unpadded, max. working volume: 2 ml, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

**Product characteristics**

Colour of print/label	green
Graduation	no
Design	unpadded
Volume of work	2 ml

**Size**

Diameter	6 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Aspiration pipette, unpadded, 2 ml, sterile,  
non-pyrogenic/endotoxin-free, non-cytotoxic, 1  
piece(s)/blister**

**Page 2**

**Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

**Purity & certification**

Product category                      no medical device | IVD  
Purity standard                        sterile, non-pyrogenic/endotoxin-free, non-cytotoxic  
Sterilisation                            Gamma irradiation  
Pyrogen-free                            yes  
Batched                                 yes

**Packaging**

Minimum order qty.                    1000  
Type of inner packaging                blister  
Piece(s) / inner box                    100  
Piece(s) / case                         1000  
Piece(s) / pallet                        54000  
Case volume                            0,0323 cbm  
Weight of product                      0,0044 kg  
Weight of case                         5,03 kg  
EAN of inner box                        4038917127202  
EAN case                                4038917043311

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**