




**Product Specification**  
**Serological pipette, plugged, 2 ml, sterile, 25**  
**piece(s)/bag**

**Page 1**



**Product description**

Order number	86.1252.025
Product description	Serological pipette, plugged, max. working volume: 2 ml, coding colour: green, material: PS, sterile, 25 piece(s)/bag

**Product characteristics**

Colour of print/label	green
Graduation	yes
Design	plugged
Volume of work	2 ml
Coding colour	green

**Size**

Diameter	6 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Serological pipette, plugged, 2 ml, sterile, 25**  
**piece(s)/bag**

**Page 2**

**Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

**Purity & certification**

Product category                      no medical device | IVD  
Purity standard                         sterile  
Sterilisation                              Gamma irradiation  
Batched                                    yes

**Packaging**

Minimum order qty.                    1000  
Type of inner packaging                bag  
Piece(s) / inner box                    25  
Piece(s) / case                         1000  
Piece(s) / pallet                        81000  
Case volume                             0,0194 cbm  
Weight of product                      0,0036 kg  
Weight of case                         4,21 kg  
EAN of inner box                        4038917127219  
EAN case                                 4038917021098

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**