




Product Specification

Serological pipette, with tip, plugged, 5 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister

Page 1



Product description

Order number	86.1253.001
Product description	Serological pipette, with tip, plugged, 5 ml, coding colour: blue, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Colour of print/label	blue
Graduation	yes
Design	with tip, plugged
Volume of work	5 ml
Coding colour	blue

Size

Diameter	8 mm
Length of product	340 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-12

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Serological pipette, with tip, plugged, 5 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation Gamma irradiation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 500
Type of smallest blister
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 27000
Depth of case 420 mm
Width of case 250 mm
Height of case 290 mm
Case volume 0.0305 cbm
Weight of product 0.008 kg
Weight of case 4.54 kg
EAN of inner box 4038917127240
EAN case 4038917021104

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-12

This document was prepared by EDP support and is valid without signature.



SARSTEDT