




Product Specification
**Serological pipette, with tip, plugged, 5 ml, sterile,
25 piece(s)/bag**

Page 1



Product description

Order number 86.1253.025
Product description Serological pipette, with tip, plugged, 5 ml, coding colour: blue, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Label/ Print with printed graduations
Colour of print/label blue/black
Graduation yes
Design with tip, plugged
Volume of work 0,6 ml - 7,5 ml
Coding colour blue

Size

Sample volume 5 ml
Diameter 8 mm
Length of product 340 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Serological pipette, with tip, plugged, 5 ml, sterile,
25 piece(s)/bag

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 25
Piece (s) / outer case 500
Piece(s) / pallet 40500
Depth of case 413 mm
Width of case 250 mm
Height of case 195 mm
Case volume 0.0201 cbm
Weight of product 0.0067 kg
Weight of case 4.03 kg
EAN of inner box 4038917127257
EAN case 4038917021111

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT