




Product Specification
Serological pipette, with tip, plugged, 10 ml,
sterile, non-pyrogenic/endotoxin-free,
non-cytotoxic, 1 piece(s)/blister

Page 1



Product description

Order number	86.1254.001
Product description	Serological pipette, with tip, plugged, 10 ml, coding colour: orange, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Colour of print/label	orange
Graduation	yes
Design	with tip, plugged
Volume of work	10 ml
Coding colour	orange

Size

Diameter	10 mm
Length of product	335 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Serological pipette, with tip, plugged, 10 ml,
sterile, non-pyrogenic/endotoxin-free,
non-cytotoxic, 1 piece(s)/blister**

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	27000
Depth of case	420 mm
Width of case	250 mm
Height of case	290 mm
Case volume	0.0305 cbm
Weight of product	0.0104 kg
Weight of case	5.74 kg
EAN of inner box	4038917127295
EAN case	4038917021128

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT