




**Product Specification**  
**Serological pipette, with tip, plugged, 10 ml,  
sterile, non-pyrogenic/endotoxin-free,  
non-cytotoxic, 1 piece(s)/blister**

**Page 1**



**Product description**

Order number	86.1254.001
Product description	Serological pipette, with tip, plugged, 10 ml, coding colour: orange, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

**Product characteristics**

Label/ Print	with printed graduations
Colour of print/label	orange/black
Graduation	yes
Design	with tip, plugged
Volume of work	1 ml - 13 ml
Coding colour	orange

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Serological pipette, with tip, plugged, 10 ml,  
sterile, non-pyrogenic/endotoxin-free,  
non-cytotoxic, 1 piece(s)/blister**

**Page 2**

**Size**

Sample volume	10 ml
Diameter	10 mm
Length of product	335 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	27000
Depth of case	420 mm
Width of case	250 mm
Height of case	290 mm
Case volume	0.0304 cbm
Weight of product	0.0104 kg
Weight of case	5.74 kg
EAN of inner box	4038917127295
EAN case	4038917021128

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-04-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**