




Product Specification
**Serological pipette, with tip, plugged, 10 ml,
sterile, 25 piece(s)/bag**

Page 1



Product description

Order number 86.1254.025
Product description Serological pipette, with tip, plugged, 10 ml, coding colour: orange, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Label/ Print with printed graduations
Colour of print/label orange/black
Graduation yes
Design with tip, plugged
Volume of work 1 ml - 13 ml
Coding colour orange

Size

Sample volume 10 ml
Diameter 10 mm
Length of product 335 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Serological pipette, with tip, plugged, 10 ml,
sterile, 25 piece(s)/bag**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 25
Piece (s) / outer case 500
Piece(s) / pallet 27000
Depth of case 413 mm
Width of case 250 mm
Height of case 270 mm
Case volume 0.0279 cbm
Weight of product 0.0091 kg
Weight of case 5.26 kg
EAN of inner box 4038917127318
EAN case 4038917021135

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT