



Product Specification
**Serological pipette, with tip, plugged, 25 ml,
sterile, 20 piece(s)/bag**

Page 1



Product description

Order number	86.1685.020
Product description	Serological pipette, with tip, plugged, 25 ml, coding colour: red, material: PS, sterile, 20 piece(s)/bag

Product characteristics

Label/ Print	with printed graduations
Colour of print/label	red/black
Graduation	yes
Design	with tip, plugged
Volume of work	3 ml - 35 ml
Coding colour	red

Size

Sample volume	25 ml
Diameter	15 mm
Length of product	345 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Serological pipette, with tip, plugged, 25 ml,
sterile, 20 piece(s)/bag**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category not regulated
CE certified no CE mark
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 200
Type of smallest bag
subpackaging
Piece(s) / inner box 20
Piece (s) / outer case 200
Piece(s) / pallet 10800
Depth of case 413 mm
Width of case 250 mm
Height of case 270 mm
Case volume 0.0279 cbm
Weight of product 0.0152 kg
Weight of case 3.61 kg
EAN of inner box 4038917133029
EAN case 4038917040785

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-05-23

This document was prepared by EDP support and is valid without signature.



SARSTEDT