




Product Specification
Demeter pipette, with tip, plugged, 1.1 ml, sterile,
25 piece(s)/bag

Page 1



Product description

Order number	86.1686.025
Product description	Demeter pipette, with tip, plugged, 1.1 ml, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Graduation	yes
Design	with tip, plugged
Volume of work	1,1 ml

Size

Diameter	4 mm
Length of product	270 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Demeter pipette, with tip, plugged, 1.1 ml, sterile,
25 piece(s)/bag

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	1000
Piece(s) / pallet	108000
Depth of case	413 mm
Width of case	250 mm
Height of case	150 mm
Case volume	0.0155 cbm
Weight of product	0.0033 kg
Weight of case	3.94 kg
EAN of inner box	4038917133043
EAN case	4038917024662

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

