



Product description

Order number	86.1686.225
Product description	Demeter pipette, without tip, plugged, 1.1 ml, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Label/ Print	with printed graduations
Colour of print/label	black
Graduation	yes
Design	without tip, plugged
Volume of work	0,5 ml - 1,1 ml

Size

Sample volume	1.1 ml
Diameter	4 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Demeter pipette, without tip, plugged, 1.1 ml,
sterile, 25 piece(s)/bag**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category not regulated
CE certified no CE mark
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest
subpackaging bag
Piece(s) / inner box 25
Piece (s) / outer case 1000
Piece(s) / pallet 108000
Depth of case 413 mm
Width of case 250 mm
Height of case 150 mm
Case volume 0.0155 cbm
Weight of product 0.0034 kg
Weight of case 4.04 kg
EAN of inner box 4038917133050
EAN case 4038917024679

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-03-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT