




**Product Specification**  
**Demeter pipette, without tip, plugged, 1.1 ml,  
sterile, 25 piece(s)/bag**

**Page 1**



**Product description**

Order number	86.1686.225
Product description	Demeter pipette, without tip, plugged, 1.1 ml, material: PS, sterile, 25 piece(s)/bag

**Product characteristics**

Graduation	yes
Design	without tip, plugged
Volume of work	1,1 ml

**Size**

Diameter	4 mm
Length of product	270 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Demeter pipette, without tip, plugged, 1.1 ml,  
sterile, 25 piece(s)/bag**

**Page 2**

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	1000
Piece(s) / pallet	108000
Depth of case	413 mm
Width of case	250 mm
Height of case	150 mm
Case volume	0.0155 cbm
Weight of product	0.0034 kg
Weight of case	4.04 kg
EAN of inner box	4038917133050
EAN case	4038917024679

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**