




**Product Specification**  
**Serological pipette, without tip, plugged, 5 ml,  
sterile, 10 piece(s)/bag**

**Page 1**



**Product description**

Order number	86.1687.010
Product description	Serological pipette, without tip, plugged, 5 ml, material: PS, sterile, 10 piece(s)/bag

**Product characteristics**

Graduation	yes
Design	without tip, plugged
Volume of work	5 ml

**Size**

Diameter	8 mm
Length of product	340 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	500
Piece(s) / pallet	40500
Depth of case	413 mm
Width of case	250 mm
Height of case	195 mm
Case volume	0.0201 cbm
Weight of product	0.0064 kg
Weight of case	3.88 kg
EAN of inner box	4038917133074
EAN case	4038917024686

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

