



Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1171.300
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 15 mm, material: LD-PE, transparent, droplet size: 35 µl - 55 µl

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3.5 ml
Width of product	15 mm
Length of product	155 mm
Drop size	35 µl - 55 µl

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,
LD-PE, transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	56000
Case volume	0.0994 cbm
Weight of product	0.0017 kg
Weight of case	9.41 kg
EAN of inner box	4038917378383
EAN case	4038917378376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

