



**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,**  
**LD-PE, transparent**

**Page 1**



**Product description**

Order number	86.1172
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, injected graduation up to 1 ml, 500 piece(s)/case

**Product characteristics**

Graduation	yes
------------	-----

**Size**

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	156 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,**  
**LD-PE, transparent**

**Page 2**

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	56000
Depth of box	261 mm
Width of box	161 mm
Height of box	231 mm
Depth of case	533 mm
Width of case	333 mm
Height of case	250 mm
Case volume	0.0444 cbm
Weight of product	0.0013 kg
Weight of case	3.47 kg
EAN of inner box	4038917085588
EAN case	4038917004220

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**