




**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 12.5 mm,**  
**LD-PE, transparent**

**Page 1**



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 86.1173.001  |
| Product description | Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, sterile, 1 piece(s)/blister |

**Product characteristics**

|            |    |
|------------|----|
| Graduation | no |
|------------|----|

**Size**

|                   |               |
|-------------------|---------------|
| Sample volume     | 3.5 ml        |
| Width of product  | 12.5 mm       |
| Length of product | 155 mm        |
| Drop size         | 30 µl - 45 µl |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 12.5 mm,**  
**LD-PE, transparent**

**Material & colours**

|                   |                                  |
|-------------------|----------------------------------|
| Product material  | Low Density Polyethylene (LD-PE) |
| Colour of product | transparent                      |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Purity standard  | sterile                 |
| Sterilisation    | Electron irradiation    |
| Batched          | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 4000          |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 4000          |
| Piece(s) / pallet             | 64000         |
| Case volume                   | 0.1309 cbm    |
| Weight of product             | 0.0023 kg     |
| Weight of case                | 11.23 kg      |
| EAN of inner box              | 4038917318341 |
| EAN case                      | 4038917318334 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

