



### **Product description**

Order number	86.1173.020
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, sterile, 20 piece(s)/bag

### **Product characteristics**

Graduation	no
------------	----

### **Size**

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	155 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 12.5 mm,**  
**LD-PE, transparent**

**Page 2**

**Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

**Purity & certification**

Product category not regulated  
CE certified no CE mark  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 4000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 4000  
Piece(s) / pallet 64000  
Depth of case 603 mm  
Width of case 438 mm  
Height of case 495 mm  
Case volume 0.1309 cbm  
Weight of product 0.0016 kg  
Weight of case 8.32 kg  
EAN of inner box 4038917318372  
EAN case 4038917318365

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-03-10

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**