




Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 12.5 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1173.300
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, 500 piece(s)/case

Product characteristics

Graduation	no
------------	----

Size

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	155 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-08-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 12.5 mm,
LD-PE, transparent

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 5000
Piece(s) / pallet 56000
Case volume 0.0994 cbm
Weight of product 0.0012 kg
Weight of case 6.89 kg
EAN of inner box 4038917378475
EAN case 4038917378468

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-08-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT