



Product description

Order number	86.1173.300
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl

Product characteristics

Graduation	no
------------	----

Size

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	155 mm
Drop size	30 µl - 45 µl

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 12.5 mm,
LD-PE, transparent

Page 2

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	56000
Case volume	0.0994 cbm
Weight of product	0.0012 kg
Weight of case	6.89 kg
EAN of inner box	4038917378475
EAN case	4038917378468

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-09-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT