




**Product Specification**  
**Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,**  
**transparent**

**Page 1**



**Product description**

Order number	86.1174.020
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 152 x 15 mm, material: LD-PE, transparent, droplet size: 40 µl - 55 µl, sterile, 20 piece(s)/bag

**Product characteristics**

Graduation	no
------------	----

**Size**

Sample volume	6 ml
Width of product	15 mm
Length of product	152 mm
Drop size	40 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,**  
**transparent**

**Page 2**

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	3200
Type of smallest subpackaging	bag
Piece(s) / inner box	400
Piece (s) / outer case	3200
Piece(s) / pallet	48000
Case volume	0.1395 cbm
Weight of product	0.0023 kg
Weight of case	8.95 kg
EAN of inner box	4038917318402
EAN case	4038917318396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**