




**Product Specification**  
**Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,**  
**transparent**

**Page 1**



**Product description**

Order number	86.1174.300
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 152 x 15 mm, material: LD-PE, transparent, droplet size: 40 µl - 55 µl, 400 piece(s)/case

**Product characteristics**

Graduation	no
------------	----

**Size**

Sample volume	6 ml
Width of product	15 mm
Length of product	152 mm
Drop size	40 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,**  
**transparent**

**Page 2**

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	400
Type of smallest subpackaging	case
Piece(s) / inner box	400
Piece (s) / outer case	4000
Piece(s) / pallet	44800
Case volume	0.0994 cbm
Weight of product	0.002 kg
Weight of case	9.07 kg
EAN of inner box	4038917378505
EAN case	4038917378499

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-11

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**