




Product Specification
Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE,
transparent

Page 1



Product description

Order number	86.1175.001
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 µl - 30 µl, sterile, 1 piece(s)/blister

Product characteristics

Graduation	no
------------	----

Size

Sample volume	6 ml
Width of product	15 mm
Length of product	146 mm
Drop size	20 µl - 30 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE,
transparent

Page 2

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	840
Type of smallest subpackaging	blister
Piece(s) / inner box	42
Piece (s) / outer case	840
Piece(s) / pallet	23520
Depth of case	533 mm
Width of case	333 mm
Height of case	250 mm
Case volume	0.0444 cbm
Weight of product	0.0016 kg
Weight of case	2.56 kg
EAN of inner box	4038917126311
EAN case	4038917020589

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT