



Product description

Order number	86.1175.020
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 µl - 30 µl, sterile, 20 piece(s)/bag

Product characteristics

Graduation	no
------------	----

Size

Sample volume	6 ml
Width of product	15 mm
Length of product	146 mm
Drop size	20 µl - 30 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE,
transparent

Page 2

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	3200
Type of smallest subpackaging	bag
Piece(s) / inner box	400
Piece (s) / outer case	3200
Piece(s) / pallet	3200
Case volume	0.1309 cbm
Weight of product	0.0022 kg
Weight of case	8.64 kg
EAN of inner box	4038917318433
EAN case	4038917318426

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT