




Product Specification
**Transfer pipette, 2 ml, (LxW): 154 x 11 mm, LD-PE,
transparent**

Page 1



Product description

Order number	86.1176.300
Product description	Transfer pipette, nominal volume: 2 ml, (LxW): 154 x 11 mm, material: LD-PE, transparent, droplet size: 35 µl - 45 µl

Product characteristics

Graduation	no
------------	----

Size

Sample volume	2 ml
Width of product	11 mm
Length of product	154 mm
Drop size	35 µl - 45 µl

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Transfer pipette, 2 ml, (LxW): 154 x 11 mm, LD-PE,
transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 5000
Piece(s) / inner box 500
Piece(s) / case 5000
Piece(s) / pallet 84000
Case volume 0.0994 cbm
Weight of product 0.0015 kg
Weight of case 8.64 kg
EAN of inner box 4038917378567
EAN case 4038917378550

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

