



### Product description

Order number	86.1176.300
Product description	Transfer pipette, nominal volume: 2 ml, (LxW): 154 x 11 mm, material: LD-PE, transparent, droplet size: 35 $\mu$ l - 45 $\mu$ l, 500 piece(s)/case

### Product characteristics

Graduation	no
------------	----

### Size

Sample volume	2 ml
Width of product	11 mm
Length of product	154 mm
Drop size	35 $\mu$ l - 45 $\mu$ l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 2 ml, (LxW): 154 x 11 mm, LD-PE,  
transparent**

**Page 2**

**Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched yes

**Packaging**

Minimum order qty. 5000  
Type of smallest case  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 5000  
Piece(s) / pallet 84000  
Case volume 0.0994 cbm  
Weight of product 0.0015 kg  
Weight of case 8.64 kg  
EAN of inner box 4038917378567  
EAN case 4038917378550

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**