




Product Specification
Transfer pipette, 2 ml, (LxW): 154 x 11 mm, LD-PE,
transparent

Page 1



Product description

Order number	86.1176
Product description	Transfer pipette, nominal volume: 2 ml, (LxW): 154 x 11 mm, material: LD-PE, transparent, droplet size: 35 µl - 45 µl, 750 piece(s)/case

Product characteristics

Graduation	no
------------	----

Size

Sample volume	2 ml
Width of product	11 mm
Length of product	154 mm
Drop size	35 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Transfer pipette, 2 ml, (LxW): 154 x 11 mm, LD-PE,
transparent**

Page 2

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	3000
Type of smallest subpackaging	case
Piece(s) / inner box	750
Piece (s) / outer case	3000
Piece(s) / pallet	84000
Depth of box	261 mm
Width of box	161 mm
Height of box	231 mm
Depth of case	533 mm
Width of case	333 mm
Height of case	250 mm
Case volume	0.0444 cbm
Weight of product	0.0012 kg
Weight of case	4.5 kg
EAN of inner box	4038917085625
EAN case	4038917004268

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT