



Product description

Order number	86.1180
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 104 x 10 mm, material: LD-PE, transparent, droplet size: 17 µl - 22 µl, 750 piece(s)/case

Product characteristics

Graduation	no
------------	----

Size

Sample volume	1 ml
Width of product	10 mm
Length of product	104 mm
Drop size	17 µl - 22 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Transfer pipette, 1 ml, (LxW): 104 x 10 mm, LD-PE,
transparent

Page 2

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	3000
Type of smallest subpackaging	case
Piece(s) / inner box	750
Piece (s) / outer case	3000
Piece(s) / pallet	252000
Depth of box	261 mm
Width of box	126 mm
Height of box	121 mm
Depth of case	533 mm
Width of case	263 mm
Height of case	140 mm
Case volume	0.0196 cbm
Weight of product	0.0005 kg
Weight of case	2.07 kg
EAN of inner box	4038917085649
EAN case	4038917004282

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT