



### **Product description**

Order number 86.1562.050  
Product description Inoculation loop, volume: 10 µl, material: PS, blue, sterile, 48 piece(s)/blister

### **Size**

Width of product 6 mm  
Length of product 200 mm

### **Material & colours**

Product material Polystyrene (PS)  
Colour of product blue

### **Purity & certification**

Product category no medical device | IVD  
Purity standard sterile  
Sterilisation Gamma irradiation  
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	blister
Piece(s) / inner box	960
Piece (s) / outer case	1920
Piece(s) / pallet	253440
Depth of box	216 mm
Width of box	161 mm
Height of box	101 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	120 mm
Case volume	0.0092 cbm
Weight of product	0.0011 kg
Weight of case	2.58 kg
EAN of inner box	4038917130882
EAN case	4038917073127

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

