



Product description

| | |
|---------------------|--|
| Order number | 86.1568.050 |
| Product description | Inoculation needle, material: PS, orange, sterile, 50 piece(s)/peel-pack |

Size

| | |
|-------------------|--------|
| Width of product | 4 mm |
| Length of product | 197 mm |

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | orange |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | sterile |
| Sterilisation | Gamma irradiation |
| Batched | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

| | |
|-------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of inner packaging | peel-pack |
| Piece(s) / inner box | 1000 |
| Piece(s) / case | 2000 |
| Piece(s) / pallet | 264000 |
| Depth of box | 216 mm |
| Width of box | 161 mm |
| Height of box | 101 mm |
| Depth of case | 338 mm |
| Width of case | 228 mm |
| Height of case | 120 mm |
| Case volume | 0,0092 cbm |
| Weight of product | 0,0012 kg |
| Weight of case | 2,79 kg |
| EAN of inner box | 4038917130974 |
| EAN case | 4038917023542 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

