



Product description

Order number	86.1569.001
Product description	Inoculation spreader, material: PS, blue, for spreading a bacterial suspension, sterile, 1 piece(s)/bag

Product characteristics

Example of use	for spreading a bacterial suspension
----------------	--------------------------------------

Size

Width of product	40 mm
Length of product	185 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	1
Piece(s) / case	500
Piece(s) / pallet	32000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0,0255 cbm
Weight of product	0,0021 kg
Weight of case	3,03 kg
EAN of inner box	4038917130981
EAN case	4038917023559

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

