



Product description

Order number	86.1569.005
Product description	Inoculation spreader, material: PS, blue, for spreading a bacterial suspension, sterile, 4 piece(s)/blister

Product characteristics

Example of use	for spreading a bacterial suspension
----------------	--------------------------------------

Size

Width of product	40 mm
Length of product	185 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	not regulated
CE certified	kein CE-Kennzeichen
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	4
Piece (s) / outer case	500
Piece(s) / pallet	33000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.002 kg
Weight of case	1.7 kg
EAN of inner box	4038917130998
EAN case	4038917023566

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

