



### **Product description**

Order number	86.1569.005
Product description	Inoculation spreader, material: PS, blue, for spreading a bacterial suspension, sterile, 4 piece(s)/peel-pack

### **Product characteristics**

Example of use	for spreading a bacterial suspension
----------------	--------------------------------------

### **Size**

Width of product	40 mm
Length of product	185 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product in accordance with MDD or IVDD	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Lot or serial number	yes

### **Packaging**

Minimum order qty.	500
Type of inner packaging	peel-pack
Piece(s) / inner box	4
Piece(s) / case	500
Piece(s) / pallet	33000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0,0182 cbm
Weight of product	0,002 kg
Weight of case	1,7 kg
EAN of inner box	4038917130998
EAN case	4038917023566

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

