



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 86.1569.005   |
| Product description | Inoculation spreader, material: PS, blue, for spreading a bacterial suspension, sterile, 4 piece(s)/blister |

### **Product characteristics**

|                |                                      |
|----------------|--------------------------------------|
| Example of use | for spreading a bacterial suspension |
|----------------|--------------------------------------|

### **Size**

|                   |        |
|-------------------|--------|
| Width of product  | 40 mm  |
| Length of product | 185 mm |

### **Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | blue             |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Purity standard  | sterile                 |
| Sterilisation    | Gamma irradiation       |
| Batched          | yes                     |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 4             |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 33000         |
| Depth of case                 | 438 mm        |
| Width of case                 | 268 mm        |
| Height of case                | 155 mm        |
| Case volume                   | 0.0182 cbm    |
| Weight of product             | 0.002 kg      |
| Weight of case                | 1.7 kg        |
| EAN of inner box              | 4038917130998 |
| EAN case                      | 4038917023566 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

