



Product Specification
DishRack, height: 370 mm, yellow, for 88 petri dishes with 92 mm Ø

Page 1



Product description

| | |
|---------------------|---|
| Order number | 93.1646.002 |
| Product description | DishRack, height: 370 mm, material: PC, yellow, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case |

Product characteristics

| | |
|------------------|----|
| Storage capacity | 88 |
|------------------|----|

Size

| | |
|-------------------|--------|
| Height of product | 370 mm |
|-------------------|--------|

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polycarbonate (PC) |
| Colour of product | yellow |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
DishRack, height: 370 mm, yellow, for 88 petri dishes with 92 mm Ø

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging case
Piece(s) / inner box 1
Piece (s) / outer case 1
Piece(s) / pallet 108
Depth of case 383 mm
Width of case 193 mm
Height of case 190 mm
Case volume 0.014 cbm
Weight of product 1.086 kg
Weight of case 1.35 kg
EAN of inner box 4038917132565
EAN case 4038917030472

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

