



Product description

Order number	93.1646.004
Product description	DishRack, height: 370 mm, material: PC, green, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

Product characteristics

Storage capacity	88
------------------	----

Size

Height of product	370 mm
-------------------	--------

Material & colours

Product material	Polycarbonate (PC)
Colour of product	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	108
Depth of case	383 mm
Width of case	193 mm
Height of case	190 mm
Case volume	0.014 cbm
Weight of product	1.086 kg
Weight of case	1.35 kg
EAN of inner box	4038917132589
EAN case	4038917030496

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

