



Product description

Order number	93.1646.080
Product description	Silicone strip, for DishRack 80, length: 370 mm, material: silicone, translucent, 1 piece(s)/bag

Product characteristics

Suitable for	DishRack 80
--------------	-------------

Size

Width of product	65 mm
Length of product	370 mm

Material & colours

Product material	silicone
Colour of product	translucent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	30
Piece(s) / pallet	3240
Depth of case	383 mm
Width of case	193 mm
Height of case	190 mm
Case volume	0.014 cbm
Weight of product	0.0579 kg
Weight of case	2 kg
EAN of inner box	4038917132596
EAN case	4038917040655

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

