



Product Specification

DishRack, height: 240 mm, blue, for 52 petri dishes with 92 mm Ø

Page 1



Product description

Order number	93.1647.003
Product description	DishRack, height: 240 mm, material: PC, blue, for 52 petri dishes with 92 mm Ø, 1 piece(s)/case

Product characteristics

Storage capacity	52
------------------	----

Size

Height of product	240 mm
-------------------	--------

Material & colours

Product material	Polycarbonate (PC)
Colour of product	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**DishRack, height: 240 mm, blue, for 52 petri
dishes with 92 mm Ø**

Page 2

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 1
Piece(s) / pallet 162
Depth of case 253 mm
Width of case 193 mm
Height of case 190 mm
Case volume 0.0093 cbm
Weight of product 0.617 kg
Weight of case 0.8 kg
EAN of inner box 4038917132626
EAN case 4038917030533

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT