



**Product Specification**  
**DishRack, height: 240 mm, transparent, for 52  
petri dishes with 92 mm Ø**

**Page 1**



**Product description**

Order number	93.1647
Product description	DishRack, height: 240 mm, material: PC, transparent, for 52 petri dishes with 92 mm Ø, 1 piece(s)/case

**Product characteristics**

Storage capacity	52
------------------	----

**Size**

Height of product	240 mm
-------------------	--------

**Material & colours**

Product material	Polycarbonate (PC)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-24

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**DishRack, height: 240 mm, transparent, for 52**  
**petri dishes with 92 mm Ø**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                    no

**Packaging**

Minimum order qty.                    1  
Type of smallest                        case  
subpackaging  
Piece(s) / inner box                    1  
Piece (s) / outer case                  1  
Piece(s) / pallet                        162  
Depth of case                            253 mm  
Width of case                            193 mm  
Height of case                           190 mm  
Case volume                            0.0093 cbm  
Weight of product                      0.617 kg  
Weight of case                          0.8 kg  
EAN of inner box                        4038917087360  
EAN case                                 4038917030502

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

