



**Product Specification**  
**Cuvette, 4.2 ml, (HxW): 45 x 12 mm, PMMA,**  
**transparent, optical sides: 2**

**Page 1**



**Product description**

Order number	67.738
Product description	Cuvette, max. working volume: 4.2 ml, material: PMMA, transparent, (HxW): 45 x 12 mm, optical sides: 2, suitable for use beginning at 300 nm, path length: 10 mm, sorted by mould cavity number, 100 piece(s)/polystyrene container

**Product characteristics**

Application	suitable for use beginning at 300 nm
Volume of work	2,5 ml - 4,2 ml
Number of optical sides	2
Path length	10 mm
Wavelength in the cuvette area	300 nm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Cuvette, 4.2 ml, (HxW): 45 x 12 mm, PMMA,**  
**transparent, optical sides: 2**

**Page 2**

**Size**

Max. volume	4.5 ml
Width of product	12 mm
Height of product	45 mm

**Material & colours**

Product material	Polymethyl methacrylate (PMMA)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	polystyrene container
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	48000
Depth of case	458 mm
Width of case	288 mm
Height of case	415 mm
Case volume	0.0547 cbm
Weight of product	0.0022 kg
Weight of case	5.74 kg
EAN of inner box	4038917083843
EAN case	4038917002721

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-10-06

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**