



Product description

| | |
|---------------------|---|
| Order number | 68.752 |
| Product description | Round cuvette, max. working volume: 3.5 ml, material: PP, height: 51 mm, 1,000 piece(s)/bag |

Product characteristics

| | |
|----------------|--------|
| Volume of work | 3,5 ml |
|----------------|--------|

Size

| | |
|-------------------|-------|
| Max. volume | 4 ml |
| Diameter | 12 mm |
| Height of product | 51 mm |

Material & colours

| | |
|------------------|--------------------|
| Product material | Polypropylene (PP) |
|------------------|--------------------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Round cuvette, 3.5 ml, height: 51 mm, PP

Page 2

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 5000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 100000 |
| Depth of case | 458 mm |
| Width of case | 408 mm |
| Height of case | 350 mm |
| Case volume | 0.0654 cbm |
| Weight of product | 0.0012 kg |
| Weight of case | 6.61 kg |
| EAN of inner box | 4038917083959 |
| EAN case | 4038917002844 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT