



Product description

Order number	82.1135.500
Product description	Petri dish, (ØxH): 35 x 10 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

Size

Max. volume	10 ml
Diameter	35 mm
Height of product	10 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Petri dish, 35 x 10 mm, transparent, with
ventilation cams**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category not regulated
CE certified no CE mark
Latex-free yes
Purity standard sterile
Sterilisation Gamma irradiation
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 500
Type of smallest
subpackaging tube bag
Piece(s) / inner box 20
Piece (s) / outer case 500
Piece(s) / pallet 30000
Depth of case 322 mm
Width of case 257 mm
Height of case 260 mm
Case volume 0.0215 cbm
Weight of product 0.0036 kg
Weight of case 2.33 kg
EAN of inner box 4038917347839
EAN case 4038917365772

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-04-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT