




**Product Specification**  
**Petri dish, 35 x 10 mm, transparent, with ventilation cams**

**Page 1**



**Product description**

Order number 82.1135.500  
Product description Petri dish, (ØxH): 35 x 10 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag

**Product characteristics**

Minimum temperature 0 °C  
Maximum temperature 80 °C  
Ventilation cams yes

**Size**

Max. volume 10 ml  
Diameter 35 mm  
Height of product 10 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-01-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Petri dish, 35 x 10 mm, transparent, with ventilation cams**

**Page 2**

**Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

**Purity & certification**

Product category                      no medical device | IVD  
Latex-free                                yes  
Purity standard                         sterile  
Sterilisation                             Gamma irradiation  
Safety data sheet                        no  
Batched                                  yes

**Packaging**

Minimum order qty.                    500  
Type of inner packaging                tube bag  
Piece(s) / inner box                    20  
Piece(s) / case                         500  
Piece(s) / pallet                        30000  
Depth of case                            322 mm  
Width of case                            257 mm  
Height of case                            245 mm  
Case volume                             0,0203 cbm  
Weight of product                       0,0036 kg  
Weight of case                          2,33 kg  
EAN of inner box                        4038917347839  
EAN case                                 4038917365772

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-01-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**