




Product Specification
**Petri dish, 60 x 15 mm, transparent, with
ventilation cams**

Page 1



Product description

Order number	82.1194.500
Product description	Petri dish, (ØxH): 60 x 15 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

Size

Max. volume	30 ml
Diameter	60 mm
Height of product	15 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Petri dish, 60 x 15 mm, transparent, with
ventilation cams**

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	not regulated
CE certified	no CE mark
Latex-free	yes
Purity standard	sterile
Sterilisation	Gamma irradiation
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	tube bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	317 mm
Width of case	317 mm
Height of case	365 mm
Case volume	0.0367 cbm
Weight of product	0.0074 kg
Weight of case	4.49 kg
EAN of inner box	4038917347891
EAN case	4038917347884

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP
support and is valid without signature.



SARSTEDT