



Product description

Order number	82.1194.500
Product description	Petri dish, (ØxH): 60 x 15 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

Size

Max. volume	30 ml
Diameter	60 mm
Height of product	15 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Petri dish, 60 x 15 mm, transparent, with ventilation cams

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Latex-free yes
Purity standard sterile
Sterilisation Gamma irradiation
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 500
Type of inner packaging tube bag
Piece(s) / inner box 20
Piece(s) / case 500
Piece(s) / pallet 12000
Depth of case 317 mm
Width of case 317 mm
Height of case 365 mm
Case volume 0,0367 cbm
Weight of product 0,0074 kg
Weight of case 4,49 kg
EAN of inner box 4038917347891
EAN case 4038917347884

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT