



Product description

Order number	82.1472.001
Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, transparent, without ventilation cams, sterile, 20 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	no

Size

Max. volume	80 ml
Diameter	92 mm
Height of product	16 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Petri dish, 92 x 16 mm, transparent, without ventilation cams

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Latex-free yes
Purity standard sterile
Sterilisation Gamma irradiation
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 480
Type of inner packaging tube bag
Piece(s) / inner box 480
Piece(s) / case 480
Piece(s) / pallet 9600
Depth of case 584 mm
Width of case 387 mm
Height of case 360 mm
Case volume 0.0814 cbm
Weight of product 0.0168 kg
Weight of case 9.19 kg
EAN of inner box 4038917129312
EAN case 4038917022118

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT