



Product Specification
Petri dish, 92 x 16 mm, green, with ventilation
cams

Page 1



Product description

| | |
|---------------------|---|
| Order number | 82.1473.080 |
| Product description | Petri dish, (ØxH): 92 x 16 mm, material: PS, green, with ventilation cams, 20 piece(s)/tube bag |

Product characteristics

| | |
|---------------------|-------|
| Minimum temperature | 0 °C |
| Maximum temperature | 80 °C |
| Ventilation cams | yes |

Size

| | |
|-------------------|-------|
| Max. volume | 80 ml |
| Diameter | 92 mm |
| Height of product | 16 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | green |

Purity & certification

| | |
|-------------------|-------------------------|
| Product category | no medical device IVD |
| Latex-free | yes |
| Safety data sheet | no |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 480 |
| Type of smallest subpackaging | tube bag |
| Piece(s) / inner box | 480 |
| Piece (s) / outer case | 480 |
| Piece(s) / pallet | 9600 |
| Depth of case | 584 mm |
| Width of case | 387 mm |
| Height of case | 360 mm |
| Case volume | 0.0814 cbm |
| Weight of product | 0.0169 kg |
| Weight of case | 9.25 kg |
| EAN of inner box | 4038917129404 |
| EAN case | 4038917022170 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

