




**Product Specification**  
**Petri dish, 100 x 100 x 20 mm, transparent,**  
**without ventilation cams**

**Page 1**



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 82.9923.422  |
| Product description | Petri dish, (LxWxH): 100 x 100 x 20 mm, material: PS, transparent, without ventilation cams, sterile, 4 piece(s)/bag |

**Product characteristics**

|                     |       |
|---------------------|-------|
| Minimum temperature | 0 °C  |
| Maximum temperature | 80 °C |
| Ventilation cams    | no    |

**Size**

|                   |        |
|-------------------|--------|
| Max. volume       | 175 ml |
| Width of product  | 100 mm |
| Height of product | 20 mm  |
| Length of product | 100 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Petri dish, 100 x 100 x 20 mm, transparent,**  
**without ventilation cams**

**Page 2**

**Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | transparent      |

**Purity & certification**

|                   |                         |
|-------------------|-------------------------|
| Product category  | no medical device   IVD |
| Latex-free        | yes                     |
| Purity standard   | sterile                 |
| Sterilisation     | Gamma irradiation       |
| Safety data sheet | no                      |
| Batched           | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 160           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 4             |
| Piece (s) / outer case        | 160           |
| Piece(s) / pallet             | 4800          |
| Depth of case                 | 608 mm        |
| Width of case                 | 324 mm        |
| Height of case                | 264 mm        |
| Case volume                   | 0.052 cbm     |
| Weight of product             | 0.0343 kg     |
| Weight of case                | 6.35 kg       |
| EAN of inner box              | 4038917154178 |
| EAN case                      | 4038917027854 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**