




Product Specification
Petri dish, 100 x 100 x 20 mm, transparent,
without ventilation cams

Page 1



Product description

Order number	82.9923.422
Product description	Petri dish, (LxWxH): 100 x 100 x 20 mm, material: PS, transparent, without ventilation cams, sterile, 4 piece(s)/bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	no

Size

Max. volume	175 ml
Width of product	100 mm
Height of product	20 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-08-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Petri dish, 100 x 100 x 20 mm, transparent,
without ventilation cams

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product in accordance with no medical device | IVD
MDD or IVDD
Latex-free yes
Purity standard sterile
Sterilisation Gamma irradiation
Safety data sheet nein
Lot or serial number yes

Packaging

Minimum order qty. 160
Type of inner packaging bag
Piece(s) / inner box 4
Piece(s) / case 160
Piece(s) / pallet 4800
Case volume 0,0577 cbm
Weight of product 0,0343 kg
Weight of case 6,35 kg
EAN of inner box 4038917154178
EAN case 4038917027854

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-08-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT