



### **Product description**

Order number	82.9923.427
Product description	Petri dish, (ØxH): 54.65 x 14.7 mm, material: PS, transparent, without ventilation cams, sterile, 10 piece(s)/bag

### **Product characteristics**

Ventilation cams	no
------------------	----

### **Size**

Diameter	54,65 mm
Height of product	14,7 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

### **Packaging**

Minimum order qty.	1000
Type of inner packaging	bag
Piece(s) / inner box	10
Piece(s) / case	1000
Piece(s) / pallet	40000
Depth of case	604 mm
Width of case	264 mm
Height of case	373 mm
Case volume	0,06 cbm
Weight of case	8,2 kg
EAN of inner box	4038917154215
EAN case	4038917039789

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

