




Product Specification
Petri dish, 68.45 x 15 mm, transparent, without ventilation cams

Page 1



Product description

Order number	82.9923.429
Product description	Petri dish, (ØxH): 68.45 x 15 mm, material: PS, transparent, without ventilation cams, sterile, 10 piece(s)/bag

Product characteristics

Ventilation cams	no
------------------	----

Size

Diameter	68.45 mm
Height of product	15 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	600
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	600
Piece(s) / pallet	28800
Depth of case	608 mm
Width of case	264 mm
Height of case	324 mm
Case volume	0.0519 cbm
Weight of case	6.7 kg
EAN of inner box	4038917154239
EAN case	4038917039802

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

