



### **Product description**

Order number	80.1301
Product description	Swabs, in the tube, no medium, (LxØ): 175 x 12 mm, stem material: PS, stem length: 133 mm, swab material: viscose, sterile, 100 piece(s)/case

### **Product characteristics**

Type of cap	push cap
Design	in the tube
Type of medium	no medium

### **Size**

Diameter	12 mm
Length of product	175 mm
Length including cap	175 mm
Length excluding cap	155 mm
Length including cap and plunger	166 mm
Length of stem	133 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	white
Swab material	viscose
Stem material	Polystyrene (PS)

### **Purity & certification**

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gamma irradiation
Safety data sheet	no
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of inner packaging	case
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	32000
Case volume	0,0367 cbm
Weight of product	0,0047 kg
Weight of case	5,94 kg
EAN of inner box	4038917086103
EAN case	4038917234184

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

